

---

# IGGPIPE INSTALLATION

Ted Toal <twtoal@ucdavis.edu>

## Table of Contents

1. If not done already, download IGGPIPE files .....	1
2. Requirements for installing IGGPIPE .....	1
3. Install IGGPIPE .....	3
4. To run IGGPIPE to generate markers .....	15
4.1. For problems and help: .....	16

This describes how to install all elements needed to run the IGGPIPE software.

## 1. If not done already, download IGGPIPE files

1. Browse to <https://github.com/BradyLab/IGGPIPE>
2. At bottom of right column on screen, click "Download ZIP" and choose a place to put it on your computer.
3. Unzip the zip file on your computer.
4. Rename the unzipped folder from "IGGPIPE-master" to just "IGGPIPE".
5. If you are a GitHub user and wish to use the IGGPIPE GitHub repository directly, instead of the above steps, you can fork the IGGPIPE repository and work on your own copy. The *master* branch is where the releases reside, each tagged with a version number. You can create a branch whose contents are the same as a tagged version with the command:

```
git checkout -b my_V1.0_branch v1.0"
```

You can then do the installation process using those files.

I will try to monitor for pull requests if you make useful changes.

## 2. Requirements for installing IGGPIPE

### 1. 64-bit Mac OSX or Linux System

This description focuses on installation on a 64-bit Mac OSX system. Installation has also been tested on 64-bit Linux and 64-bit Windows, and we describe installation differences, but command line examples here are for Mac OSX, and Linux and Windows users are expected to be savvy enough to translate these to their system. IGGPIPE will not work on 32-bit Windows systems.

### 2. The computer must have enough memory and speed

Before starting the installation, consider whether the computer on which the software is to be installed has enough memory for the problem at hand. I have been running IGGPIPE using 16GB of RAM, and I'd recommend no less than that amount. For larger genomes or more genomes, you will possibly need more. My computer has a 2.4 GHz processor, which has been adequate. You should monitor memory usage while the pipeline is running, and if it becomes excessive for your computer, you know you will need more memory to run those genomes.

### 3. Work from the command line

Most work here is done from the command line, by opening the Terminal application. You should be familiar with some basic command line commands such as *cd*, *ls*, *cp*, *rm*, *mv*, *sudo*, *less* or *more*, *head* and *tail*, *man*, and *mkdir*. If you do not know how to use the command line or don't know these basic commands, take time to learn a bit. A suggested tutorial (on the long side, but spend as little or as much time as you need to feel you can get started here) from UC Davis is at:

[http://korflab.ucdavis.edu/Unix\\_and\\_Perl/current.html](http://korflab.ucdavis.edu/Unix_and_Perl/current.html)

### 4. Work in the IGGPIPE main directory unless otherwise instructed

While working from the command line to install IGGPIPE, most of the time you will be in the IGGPIPE main directory, unless instructed otherwise. If you unzipped the IGGPIPE zip file in your Documents folder, you might change into the IGGPIPE directory with this command:

```
cd ~/Documents/IGGPIPE
```

### 5. Know how to use a plain text editor and have one available

You must have a plain text editor you know how to use. If nothing else, the simple editor *nano*, which comes with Mac and Linux systems, will work. On a Mac you can also use the *TextEdit* program (use Plain Text format). For Mac, the free open-source TextWrangler program is strongly recommended, and its more capable version called BBEdit is totally worth the \$50 price for anyone who regularly edits text files. TextWrangler is available from the Apple App Store (Applications, App Store) or from:

<http://www.barebones.com/products/textwrangler>

### 6. The following programs/applications will be installed if they aren't already:

- a. Xcode (V6.4 works) (Apple App Store; *make* and *g++*)
- b. Jellyfish V2 or later (V2.2.3 works) (<http://www.genome.umd.edu/jellyfish.html>)
- c. Perl V5 or later (V5.16.0 works) (<https://www.perl.org/get.html>)
- d. R V3 or later (V3.2.1 works) (<https://www.r-project.org>)
- e. Primer3 (V2.3.6 works) (<http://sourceforge.net/projects/primer3>)
- f. NCBI e-PCR (V2.3.12-1 works) (<http://www.ncbi.nlm.nih.gov/tools/epcr>)
- g. MUSCLE (V3.8.31 works) (<http://www.drive5.com/muscle/downloads.htm>)

### 7. Two of the above will require compilation and build on your computer using make:

- a. Jellyfish
  - b. e-PCR
8. **One C++ program provided with IGGPIPE will be compiled and built:**
- a. findMers

## 3. Install IGGPIPE

**If you want a Quick Start with minimal reading, just look for the command lines below and execute them, with occasional skimming of the text just above/below the commands.**

### 1. Change to IGGPIPE directory

All work from the command line while installing IGGPIPE is generally done with the current directory set to the IGGPIPE directory where you download the IGGPIPE files:

```
cd whatever/IGGPIPE
```

In Windows under Cygwin, you must use the following command after changing to the IGGPIPE directory:

```
cygpath -am
```

### 2. Copy allPathParameters.template file to new file allPathParameters.ours

A variety of applications must be installed. Before installing them, you should prepare a text file for editing. There is a text file in the IGGPIPE main directory that contains tool path settings for running IGGPIPE: allPathParameters.template.

The file provides settings of paths where applications have been installed, and related settings. To make your own version of this file containing your own application paths, copy the file to a new filename, replacing ".template" with ".ours":

```
cp allPathParameters.template allPathParameters.ours
```

### 3. Open the allPathParameters.ours file in a plain text editor

Open the new allPathParameters.ours file created above in your plain text editor for editing. If you are in a hurry, you don't need to go through the whole file, but simply need to set the parameters shown at the start of the file, up to the comment that indicates you are at the end of the quick start section.

For example, if your text editor is nano, you might use this command line to open your editor to edit the template file:

```
nano allPathParameters.ours
```

### 4. Set WD to your IGGPIPE directory in the allPathParameters.ours file

Find the WD parameter in the allPathParameters.ours file, which looks like:

```
WD := $(BRADYLAB)/Genomes/kmers/IGGPIPE
```

Change the assigned value to the path of your IGGPIPE directory (where you unzipped the IGGPIPE files). For example, maybe it would look like this:

```
WD := /Users/johndoe/Documents/IGGPIPE
```

## 5. Install *make* and *g++* (Xcode in Mac OSX, Cygwin in Windows)

IGGPIPE makes use of a utility called *make*, and also, some of the applications used by IGGPIPE are distributed as source code that must be compiled and built into a runnable application on the user's computer, which requires a C++ compiler (*g++* utility). On Linux, these utilities are already installed so you can skip this step. On Mac OSX, the Apple Developer Toolkit named Xcode provides these utilities, and it is available free from the Apple App Store (Applications, App Store). On Windows, these utilities are provided by the Cygwin tools, available from:

```
https://www.cygwin.com
```

For Mac OSX, if you don't have Xcode installed already, run the App Store application, search for "Xcode", and double-click the *Install* button to install it, and even if you do have it installed, make sure you are updated with the latest version. I used version 6.4, although later versions should work fine. Installation takes quite a long time, during which it appears nothing is happening. When it is finished, you can verify that it was installed successfully by finding the Xcode application icon in Applications and running it. It may then display a box requesting your computer administrator password so it can install additional components. Then, close the Xcode application and go to the command line and enter the following command, which checks to see if the command line tools such as *make* and *g++* are installed, and if not, installs them:

```
xcode-select --install
```

To verify they are installed, you can enter this command:

```
g++
```

and you should see the error message "clang: error: no input files".

## 6. Install Jellyfish

Jellyfish is a free open-source bioinformatics application that searches FASTA sequence files for k-mers of a specified size and writes them to a file. IGGPIPE uses Jellyfish to extract unique (occurring once) k-mers from the genome sequences being used. You may already have Jellyfish installed, and if so you want to check its version number. Here is the command to see if it is installed and check the version:

```
jellyfish --version
```

Assuming you do not have it installed, or you have a version older than 2.2.3, you must do an installation. You can find the Jellyfish at:

```
http://www.genome.umd.edu/jellyfish.html
```

I chose the "latest source and binaries" link, then downloaded the .tar.gz file. I double-clicked this file in Finder, in the Downloads folder, and it unpacked to produce a jellyfish folder. I moved this folder to a directory I made named *src* under my user root directory:

```
cd ~
```

```
pwd
mkdir src
cp Downloads/jellyfish-2.2.3 src
```

This version of IGGPIPE was tested with Jellyfish version 2.2.3. Newer versions should work as well. *Older versions will not work, because Jellyfish changed its output file names. They used to end with "\_0" but no longer do!*

Now the jellyfish program must be compiled and built into an application, and installed on your computer. I used these commands:

```
cd ~/src/jellyfish-2.2.3
./configure
make
```

The commands worked without error on OSX and Linux, but compile errors ("impossible constraint" errors) occurred on Windows. We fixed these by editing file `file_header.hpp` and adding the following lines *after* `#include <jellyfish/rectangular_binary_matrix.hpp>`:

```
namespace std {
    #include <sstream>

    template <typename T>
    std::string to_string(T value)
    {
        //create an output string stream
        std::ostringstream os ;

        //throw the value into the string stream
        os << value ;

        //convert the string stream into a string and return
        return os.str() ;
    }

    template std::string to_string<long long int>(long long int); // instantiate with
}
```

That allowed the *make* to succeed under Windows too.

Then, to install the code, use this command:

```
sudo make install
```

The *sudo* command prompts for a password, and I entered my computer's administrator password. The above command places the Jellyfish program in the PATH variable so that the program can be run with the command *jellyfish*. (Sometimes it is necessary to log out and back in for a PATH change to take effect). Rerun the "*jellyfish --version*" command to verify that Jellyfish is installed.

The parameter `CMD_JELLYFISH` in the `allPathParameters.ours` file is already set to *jellyfish*, which is the command needed to run the Jellyfish program. You shouldn't need to change it.

Also, you shouldn't need to change the parameter `JELLYFISH_HASH_SIZE`. The value that is set for it already will usually work fine. However, if you are using a computer with lots of memory, you may want to change the value to take advantage of that. It can be especially helpful if you are working with k-mer sizes or genome sizes that produce lots more than 25 million k-mers.

## 7. Install Perl

Perl is a programming language used by IGGPIPE. Using it requires a Perl interpreter application on your computer. The Mac OSX system comes with a Perl interpreter already installed, and this should be sufficient. This version of IGGPIPE was tested with Perl version 5.16.0, although later versions, and earlier V5 versions, will probably be fine. You can find out if you already have Perl installed and what its version is with this command:

```
perl --version
```

If you do not have Perl installed or if the version is older than V5, you must install it, so look for it here:

<https://www.perl.org/get.html>

Explicit installation instructions are not given here. Follow the instructions provided in the downloaded installation package, then re-run the "perl --version" command to verify that it is installed. Sometimes it may be necessary to log out and log back in so that the Perl location can be added to the PATH, before the command will work.

The parameter CMD\_PERL in the allPathParameters.ours file is already set to *perl*, which is the command needed to run the Perl program. You shouldn't need to change it.

## 8. Install R

R is a programming language used by IGGPIPE. Using it requires that the R programming environment be installed on your computer. This version of IGGPIPE was tested with R version 3.2.1, although later versions, and earlier V3 versions, will probably be fine. You can find out if you already have R installed and what its version is with this command, which invokes the command line version of the R interpreter:

```
Rscript --version
```

If you do not have R installed, or have it installed but want to update to a newer version number, look for it here:

<https://www.r-project.org>

Explicit installation instructions are not given here. Follow the instructions provided in the downloaded installation package, then re-run the "Rscript --version" command to verify that it is installed.

IGGPIPE does not use any extra R packages.

The parameter CMD\_RSCRIPT in the allPathParameters.ours file is already set to *Rscript*, which is the command needed to run the Rscript program. You shouldn't need to change it.

## 9. Install Primer3

Primer3 is a classic bioinformatics application that generates primers from sequence data. It is used by IGGPIPE to generate primers for candidate IGG markers, so it must be installed on your computer. This version of IGGPIPE was tested with Primer3 version 2.3.6, although later versions

and earlier V2 versions will probably be fine. You probably know if you already have Primer3 installed. If you don't know that you do, then you should install it. Look for it here:

<http://sourceforge.net/projects/primer3>

It comes pre-built for OSX and Windows but may need to be compiled for Linux. Make sure you download the correct version (primer3, not primer3plus). Put the downloaded directory wherever you want on your computer. The file named primer3\_core (primer3\_core.exe on Windows) in the root directory of the downloaded package is the executable program file.

Now assign the parameter "CMD\_PRIMER3CORE" in the allPathParameters.ours file, for example:

```
CMD_PRIMER3CORE := ~/Documents/primer3-2.3.6/primer3_core
```

## 10. Install e-PCR

e-PCR is an "electronic PCR" application from NCBI that uses primers and sequence data to do an *in silico* PCR amplification. It is used by IGGPIPE to test primers of candidate IGG markers to see if they generate unique amplicons of the expected length, so it must be installed on your computer. This version of IGGPIPE was tested with e-PCR version 2.3.12-1 (-V option displays version 2.3.12, but downloaded file was 2.3.12-1), although later versions will probably be fine.

To install e-PCR, look for it here:

<http://www.ncbi.nlm.nih.gov/tools/epcr>

The download link uses FTP protocol. Log in as user GUEST with no password. Look for the latest .zip version or tar.gz, copy the file or folder to your computer, and unzip it. Put the unzipped directory wherever you want on your computer.

In some cases, a binary distribution might be available, so once downloaded, you should be able to run e-PCR without further ado. At the time we downloaded version 2.3.12-1, it was only available as source code and it was necessary to run *make* to compile and build the program.

Version 2.3.12-1 had two problems with it that required editing of the source code in order for the *make* operation to complete successfully under OSX. Perhaps these problems will have been fixed in the version you download (or perhaps a binary version will be available at the time you download). Test by trying to build e-PCR. Refer to the file *BUILD.html* in the e-PCR source directory for instructions on compiling the source. For Mac OSX, the source was compiled by changing into the directory where the files were unzipped and entering the following command:

```
cd e-PCR-2.3.12-1
make LF64LDFLAGS= LF64CCFLAGS=-DNATIVE_LARGEFILES COMMON_CC_FLAGS=-w
```

If the *make* completes without error, there will be a file named "e-PCR" in the directory, and if you run it, it will display a page full of usage info:

```
e-PCR      (Run e-PCR to see if it works)
```

If you get errors from the *make* like I did, here are the changes I made that allowed the *make* to succeed:

- a. Edit file mmap.cpp and remove `"/"` from the start of the line that reads `"/#include <sstream>"`

- b. Edit file `minilcs.hpp` and insert the following two lines after the line that reads `"#include <cstring>":`

```
#include <cstdlib>
#include <sstream>
```

Now try the *make* command again, followed by running "e-PCR":

```
make LF64LDFLAGS= LF64CCFLAGS=-DNATIVE_LARGEFILES COMMON_CC_FLAGS=-w
e-PCR      (Run e-PCR)
```

The *make* should succeed and e-PCR should display its usage information, meaning you are good to go.

Now assign the parameter "CMD\_EPCR" in the `allPathParameters.ours` file, for example:

```
CMD_EPCR := ~/Documents/e-PCR-2.3.12-1/e-PCR
```

## 11. Install MUSCLE

MUSCLE is an open-domain multiple sequence aligner. It is used by IGGPIPE only if you choose to search markers or LCRs for Indels by using the *make Indels* command, so if you don't do that you can skip this step, although you may as well install it. This version of IGGPIPE was tested with MUSCLE version v.8.31, although later versions will probably be fine. To install MUSCLE, look for it here:

<http://www.drive5.com/muscle/downloads.htm>

The executable images are already built, so choose the correct download for your system and download the file, putting it wherever you want on your computer, such as a bin folder.

Now assign the parameter "CMD\_ALIGNER" in the `allPathParameters.ours` file, for example:

```
CMD_ALIGNER := ~/bin/muscle3.8.31_i86darwin64
```

## 12. Build findMers

`findMers` is a C++ program that is part of IGGPIPE. It takes as input a file full of k-mers and a genome FASTA file, and produces as output a file of the k-mers with their genomic position included as additional data columns in the file. It can also locate all contigs in the genome FASTA file and output a file that lists the starting position and length of each contig. IGGPIPE uses both of these functions of `findMers` to generate a list of common unique k-mers to be analyzed for LCRs (locally conserved regions). The `findMers` program must be compiled and built using *make*. Its source files are located in the `code/cpp/findMers` directory. Change into that directory and enter the command *make*:

```
cd code/cpp/findMers
make
findMers
cd ../../..
```

The *make* should compile the C++ files in the `findMers` folder. It should complete without error, and there will be a file named "findMers" in the directory, and when that file is run with the



*findMers* command shown above following *make*, it will display a page of usage information. The path to "findMers" is already set correctly in the `allPathParameters.ours` file.

### 13. Test trashing and choose deletion method

IGGPIPE uses *make* to run data through its pipeline. A command can be given to cause *make* to delete files that it has generated by running the pipeline. There are two different ways it can delete files: it can actually delete them, or it can move them to a trash folder where they can be found and undeleted if necessary. A script file (`code/shell/trash.sh`) is provided to move files to the Mac OSX trash folder, but for linux or Windows, you must either modify that script file so that it will work with your operating system, or choose the other method that simply deletes files.

You must choose which of these methods you want. Since the trash folder method is more useful and flexible, it is the default method, but again, on Linux or Windows you will need to change it or modify `trash.sh` to work properly.

You select the method by setting the `allPathParameters.ours` parameter `CMD_DELETE_WHEN_CLEANING` to either `$(CMD_DELETE)` or `$(CMD_TRASH)`. You should make sure it is set the way you want. Also, you should test the shell script that moves files to the trash, to make sure it works. To do this, use these commands:

```
cp help.txt junk.txt
$SHELL code/shell/trash.sh junk.txt
```

Now look in the trash can to see if file "junk.txt" is there. If this doesn't work, you should set the `$(CMD_DELETE)` method as the delete method:

```
CMD_DELETE_WHEN_CLEANING := $(CMD_DELETE)
```

### 14. Copy `primer3settings.default.txt`

Primer3 uses a settings file to control many of the settings it uses to generate primers. Several sample settings files come with Primer3, in its root directory. One of these, **`primer3web_v4_0_0_default_settings.txt`**, was copied and modified for use with IGGPIPE. The file is named **`primer3settings.default.txt`**, in the main IGGPIPE directory. The following required changes were made to it:

- a. `P3_FILE_ID` was set to a descriptive settings title.
- b. `PRIMER_EXPLAIN_FLAG` was changed from 1 to 0.
- c. `PRIMER_PRODUCT_SIZE_RANGE` was set to a simplified 36-300 (primers are designed with most intervening DNA sequence removed)
- d. `PRIMER_NUM_RETURN` was changed from 5 to 1.
- e. `PRIMER_GC_CLAMP` was changed from 0 to 1 (optional but recommended).

You need to copy the default settings file to a new file that can be edited by you, should you want to change Primer3 settings for your needs while keeping a pristine copy in the original `primer3settings.default.txt` file. Copy it to this file name:

---

```
cp primer3settings.default.txt primer3settings.txt
```

This file copy is all you need to do, IGGPIPE will work with this version, and this is the required version for running the test of IGGPIPE.

The RUN instructions for IGGPIPE indicate that primer3settings.txt should be edited if you want to change primer settings for your needs. However, whenever you want to run the test of IGGPIPE as shown below, you should re-do the above copy to use the pristine file for testing.

### 15.Enable Access to FileMerge (optional and Mac only)

Parameter settings files (allParameters.\* and allPathParameters.\*) and Primer3 settings files (primer3settings.txt) can be edited by the user. You might at some time wish to see what changes were made to a file by comparing it to another similar file. The *diff* command can be used on the command line to do this. Another program, available on Mac OSX, is *FileMerge*, a great file comparison and merging tool that comes with Xcode. It is initially hidden within Xcode, but you can put it in your dock to make it more easily accessible.

To run FileMerge, start Xcode, then on the menu choose Xcode, Open Developer Tool, FileMerge. When it opens up, find its icon on the dock and set it to stay put in the dock, then you can close Xcode and in the future get to it directly from the dock.

When you run FileMerge, it prompts for two or three or four file names. To see an example of use, enter the first two file names, "left" and "right", setting "left" to allParameters.template and "right" to allParameters.test, then click "Compare". You will see a comparison of the two files, with the differences clearly shown. If you wanted to incorporate changes from one of these files into the other, you can do this easily by using the up/down arrow keys to go through the differences one by one, and use the left/right arrow keys to select whether you want the left or right side file text in the output, and you can also click in the box on the bottom that shows the merged text and edit it; when finished you can save the merged text to a new file or overwrite one of the two compared files, using File, Save Merge. Since we don't want to merge these files, exit FileMerge without saving anything.

### 16.Run IGGPIPE using the test parameters in allParameters.test and check for success

Everything is now ready to run the IGGPIPE pipeline. Data for testing it is provided in the testFASTA folder. This consists of two FASTA files that are truncated versions of the *S. lycopersicum* (tomato) and *S. pennellii* genomes, with only two chromosomes (1 and 2) and only about 14 Mbp for each one. The parameter file allParameters.test has parameters set for using these FASTA files and doing the test. It is more-or-less a copy of the allParameters.template file, modified for testing IGGPIPE.

To test IGGPIPE, from the command line in the IGGPIPE main directory, enter this command:

```
make PARAMS=allParameters.test ALL | tee logFiles/makeLog.test.txt
```

If all goes well, the pipeline will run quickly, and after four or five minutes, it should finish with the message **ALL files are up to date**.

The *tee* command routes the piped log output from *make* to the console and to the file logFiles/makeLog.test.txt. You can examine this file after the run to see what specifically happened at each step, for example with this command:

```
more logFiles/makeLog.test.txt
```

Note that the output includes timestamps telling how long each step took to run.

If the pipeline fails, an error message of some kind is displayed, and *make* stops. (There is a problem with Windows, where sometimes *make* does not stop on an error, but keeps going. We have not found a way around this. If this happens to you, you will need to go back through the output to look for errors.) If an error occurs, proceed to the next step, troubleshooting.

If no error occurs, there should be several files in the output folder "outTestHP11", including files starting with these prefixes and suffixes (shown in the order that they are produced by the pipeline):

- a. LCRs\_\*.tsv
- b. BadKmers\_\*.tsv
- c. IndelGroupsOverlapping\_\*.tsv
- d. IndelGroupsNonoverlapping\_\*.tsv
- e. NonvalidatedMarkers\_\*.tsv
- f. MarkerErrors\_\*.tsv
- g. MarkersOverlapping\_\*.tsv
- h. MarkersNonoverlapping\_\*.tsv
- i. MarkerCounts\_\*.pdf
- j. MarkerDensity\_\*.png

The MarkersOverlapping\_ and MarkersNonoverlapping\_ files are the final output files containing the markers.

The .pdf and .png files should be examined to see how they depict marker counts and densities.

The tables at the end of the RUN document describe the columns in these tab-separated data files.

To make sure the pipeline ran correctly, compare the MarkersOverlapping\_ file to the expected result, which is in subdirectory outTestHP11/goodTest:

```
diff outTestHP11/MarkersOverlapping_*.tsv outTestHP11/goodTest/
```

This command should not produce any output, indicating the two files are identical. If it produces output indicating non-identity of the files, you have a problem, so proceed to the next step, troubleshooting.

## 17. Troubleshooting

A common problem is with file paths. Pay close attention to error messages at the end before *make* stops. Recheck file paths if messages indicate a file could not be found. Note that with Windows,

which uses "\" rather than "/" to separate directories in file paths, we found that we could use "/" in all the paths in the allParameters.ours file and allParameters.test file and it worked fine; we did not have to use "\" anywhere.

Windows gave the most problems, and the most common problem with Windows was in text file line endings, which under Windows can be either "DOS" or "Unix" line endings. IGGPIPE produces files with Unix line endings exclusively, but it generally tolerates input files with either type of line ending. Most tools and programs you might use to examine the files will also tolerate either type of line ending, but occasionally, a program requires DOS line endings. Be aware of this situation during troubleshooting, and consider whether the observed problem might be one with line endings.

Another problem can be program versions. If you use an older or newer version of a program than what we used, the pipeline might fail, depending on what the changes are, or it might produce different output. Look carefully at version numbers and check to see if the output differs for any program that has a different version number than what we used.

If IGGPIPE produces a different marker output file than expected, as indicated by output being produced by the *diff* command shown in the preceding step, you should do difference testing on other output files. Each of the files whose prefixes and suffixes are listed in the previous step have a "good" version of the file containing the expected results, in folder outTestHP11/goodTest. Each of those files can be compared to the output IGGPIPE produced when you ran it using a *diff* command to see which ones are good. No output means the files match. A shell script named *diffKeyFiles.sh* is provided that runs *diff* on each of these files. To use it with the allParameters.test output files:

```
source code/shell/diffKeyFiles.sh outTestHP11 goodTest
```

It will show only a single line of output for each file, saying it diffed the file, if the files match. If they don't match, you will get a lot of output from the mismatches. A single file can be diffed with this command, for example to diff the LCRs\_ file:

```
diff outTestHP11/LCRs_*.tsv outTestHP11/goodTest/
```

If the final output file does not match, but one or more output files do match (starting with the first file listed in the previous step), then you can tell which step produced an incorrect result based on which file in the list is the first one that is incorrect. The following *make* steps produce the following output files (italicized output files are those available in the goodTest subdirectory for comparison to your files):

<b><i>make</i> Command</b>	<b>Produces output file(s)</b>
a. <i>make</i> PARAMS=myFilename getSeqInfo GENOME=ALL	GenomeData/*.idlens
b. <i>make</i> PARAMS=myFilename getContigFile GENOME=ALL	GenomeData/*.contigs
c. <i>make</i> PARAMS=myFilename getKmers GENOME=ALL	Kmers/Kmers_*.kmers
d. <i>make</i> PARAMS=myFilename kmerStats GENOME=ALL	Kmers/Kmers_*.stats

<b><i>make</i> Command</b>	<b>Produces output file(s)</b>
e. <code>make PARAMS=myFilename sortKmers GENOME=ALL</code>	<code>Kmers/Kmers_*.sorted</code>
f. <code>make PARAMS=myFilename kmerIsect</code>	<code>Kmers/isect.kmers</code>
g. <code>make PARAMS=myFilename getGenomicPos GENOME=ALL</code>	<code>Kmers/Kmers_*.isect</code>
h. <code>make PARAMS=myFilename mergeKmers GENOME=ALL</code>	<code>Kmers/Kmers_*.merge</code>
i. <code>make PARAMS=myFilename getCommonUniqueKmers</code>	<code>Kmers/common.unique.kmers</code>
j. <code>make PARAMS=myFilename findLCRs</code>	<code>LCRs_*.tsv, BadKmers_*.tsv</code>
k. <code>make PARAMS=myFilename findIndelGroups</code>	<code>IndelGroupsOverlapping_*.tsv, IndelGroupsNonoverlapping_*.tsv</code>
l. <code>make PARAMS=myFilename getDNaseqs GENOME=ALL</code>	<code>DNaseqs_*.dnaseqs</code>
m. <code>make PARAMS=myFilename findPrimers</code>	<code>NonvalidatedMarkers_*.tsv</code>
n. <code>make PARAMS=myFilename ePCRtesting GENOME=ALL</code>	<code>MarkerErrors_*.tsv</code>
o. <code>make PARAMS=myFilename removeBadMarkers</code>	<code>MarkersOverlapping_*.tsv, MarkersNonoverlapping_*.tsv</code>
p. <code>make PARAMS=myFilename plotMarkers</code>	<code>MarkerCounts_*.pdf, MarkerDensity_*.png</code>

The different results might be due to running a software package of a different version than what we used for testing. For example, a different version of e-PCR might cause a mismatch starting at file *MarkerErrors\_\*.tsv*. Another possibility to be aware of is that the files may in fact be identical except one might have Unix line ends and the other might have DOS line ends.

## 18. Run make Indels to align markers and find Indels

An R program that is NOT run as part of the pipeline when the *make ... ALL* target is built, but which can be run using *make ... Indels*, is able to read a file of LCRs, non-overlapping IndelGroups, or non-overlapping Markers, extract the DNA sequences from the genomes in each LCR or Marker region and align them, then locate all Indels in the aligned sequences and write their positions to a file. The input file to be used is selected with parameter `PATH_INDELS_INPUT_FILE` in the `allParameters.*` file. The program to find Indels is called `alignAndGetIndels.R`. Run it as follows:

```
make PARAMS=allParameters.test Indels
```

Check that the output file exists with:

```
ls outTestHP11/Markers*.indels.tsv
```

This should list the file *outTestHP11/*

*MarkersNonoverlapping\_K11k2L100D10\_2000A100\_2000d10\_100N2F0X20V3000W8M3G1.indels.tsv*

You can examine it with Excel or a text editor to see the Indel data it contains.

## 19. Run make plotIndels to plot Indel information

Another R program that is NOT run as part of the pipeline when the *make ... ALL* target is built, but which can be run using *make ... plotIndels*, reads the Indels file produced by *make ... Indels* and plots information from it in a pdf file. The program is called *plotIndels.R*. Run it as follows:

```
make PARAMS=allParameters.test plotIndels
```

Check that the output file exists with:

```
ls outTestHP11/Markers*.indels.pdf
```

This should list the file *outTestHP11/*

*MarkersNonoverlapping\_K11k2L100D10\_2000A100\_2000d10\_100N2F0X20V3000W8M3G1.indels.pdf*

You might want to open it and look at the plots.

## 20. Run dotplot.R to make a dot plot

The *LCRs\_* file contains a list of common unique k-mers assigned to locally conserved regions (LCRs), and it can be used to make a dotplot depicting alignment of the two genomes. The R program *dotplot.R* is provided to do this. It is driven by a parameter file, a sample of which has been provided, *dotplot.template*, that is set for using the test data produced by running IGGPIPE with *allParameters.test*. Run *dotplot.R* as follows:

```
Rscript code/R/dotplot.R dotplot.template
```

Check that the output file exists with:

```
ls outTestHP11/LCRs_*.dotplot.png
```

This should list the file *outTestHP11/LCRs\_K11k2L100D10\_2000.dotplot.png*, an image file. You may want to examine it (e.g. in the OSX Preview app) to see the dot plot.

There are other sample parameter files in subdirectory *dotplot*, although the parameter file is fairly straightforward and you probably don't need other examples to work from.

## 21. Run annotateFile.R to make new files containing annotated marker data in different formats

A common need is to add additional annotation information the table of markers. For example, you might be working with an introgression line population and wish to annotate each marker with the names of the lines whose introgressions that marker lies within, along with the marker position relative to the introgression. Or, you might want to annotate each marker with the ID of the nearest gene and its distance away. You may also want to change file format, from .tsv (tab-separated) to .gff3 or .gtf for adding the markers to a browser track. All this can be done with the R program *annotateFile.R* that is provided with IGGPIPE. It is driven by a parameter file, a sample of which has been provided, *annotate.template*, that is set for using the test data produced by running IGGPIPE with *allParameters.test*, along with additional annotation test data in folder *code/R/test\_GFFfuncsAndMergeData*. Run *annotateFile.R* as follows:

```
Rscript code/R/annotateFile.R annotate.template
```

Check that the output file exists with:

```
ls outTestHP11/MarkersAnnotated.*
```

This should list file *MarkersAnnotated.test.tsv* in the outTestHP11 folder. You can examine this file with a text editor or Excel to see the new column that was added compared to the input file *MarkersOverlapping\_K11k2L100D10\_2000A100\_2000d10\_100N2F0X20V3000W8M3G1.tsv*.

There are other sample parameter files in subdirectory *annotate* which produce other types of files or do other types of file data manipulation. The parameters can be challenging to set properly, especially when merging data from a separate file, so these sample files can be helpful. Also, when .gff3 files are used, they must conform well to the expected GFF3 format else an error is likely to occur.

## 22. Edit primer3settings.txt (optional)

After finishing installation, and prior to any run of IGGPIPE, you may want to edit primer3settings.txt file and make any changes that are important for your needs. For example, you might change the parameters that determine the acceptable *range of primer Tm values*. If you have several different setting values you use, you will probably want to keep a directory of different primer3settings.txt files and copy the needed one prior to each run of IGGPIPE.

The Primer3 user manual ([http://primer3.sourceforge.net/primer3\\_manual.htm](http://primer3.sourceforge.net/primer3_manual.htm)) describes all the parameters.

An explanation of the sequence data IGGPIPE gives Primer3 in order to generate primers will be helpful, particularly in understanding the setting of the parameter PRIMER\_PRODUCT\_SIZE\_RANGE. Since IGGPIPE is making primers to be used in different genomes with different sequences and sequence lengths between the two primer sites, it cannot use the typical method of giving Primer3 the entire sequence between the two primer sites. Instead, IGGPIPE gives Primer3 the concatenation of two short sequences, one around each of the two k-mers that define and anchor the candidate IGG marker. Each sequence is equal to K plus twice EXTENSION\_LEN in length. Both K (the k-mer length) and EXTENSION\_LEN (the number of bases to add on each side of the k-mer) are defined in allParameters.template. Thus, the sequence that Primer3 uses for designing the primers is equal to 2K + 4\*EXTENSION\_LEN in length. IGGPIPE also gives Primer3 a value for its parameter SEQUENCE\_PRIMER\_PAIR\_OK\_REGION\_LIST. This tells Primer3 to design one primer in the left half of the sequence and one primer in the right half. Thus, the primer product size will appear to Primer3 to be much smaller than the actual amplicon size will be, which is why PRIMER\_PRODUCT\_SIZE\_RANGE can be set to a smaller value than the amplicon sizes.

Although Primer3 is a stable program and unlikely to change a lot, if new versions of Primer3 add parameters, you might want to incorporate them into primer3settings.txt. You will see new parameters if you compare primer3settings.txt to Primer3's file primer3web\_v4\_0\_0\_default\_settings.txt (for example by using *diff* or *FileMerge*).

**That completes the installation of IGGPIPE.**

# 4. To run IGGPIPE to generate markers

- Find file RUN.pdf or RUN.html in the IGGPIPE folder on your computer and open either one and follow the instructions.

## 4.1. For problems and help:

- Post an issue on GitHub under BradyLab/IGGPIPE repository
- Contact me, Ted Toal, [twtoal@ucdavis.edu](mailto:twtoal@ucdavis.edu) [<mailto:twtoal@ucdavis.edu>]