
RUNNING IGGPIPE TO GENERATE IGG MARKERS

Ted Toal <twtoal@ucdavis.edu>

Table of Contents

1. If not done already, download IGGPIPE files	1
2. If not done already, install IGGPIPE	1
3. Requirements for running IGGPIPE	1
4. Running IGGPIPE	2
5. Post-processing tools	12
5.1. For problems and help:	14
6. Tables	14

This describes how to run the IGGPIPE software to analyze two or more genome sequences and produce a file of candidate IGG marker primer pairs for amplifying length-polymorphic regions of those genomes.

1. If not done already, download IGGPIPE files

1. Browse to <https://github.com/BradyLab/IGGPIPE>
2. At bottom of right column on screen, click "Download ZIP" and choose a place to put it on your computer.
3. Unzip the zip file on your computer.
4. Rename the unzipped folder from "IGGPIPE-master" to just "IGGPIPE".

2. If not done already, install IGGPIPE

- Look inside the downloaded IGGPIPE folder on your computer for file INSTALL.pdf or INSTALL.html and open either one and follow the instructions.

3. Requirements for running IGGPIPE

1. Obtain genome FASTA files

IGGPIPE uses as input two or more FASTA files of complete genome sequences (at least one of which is assembled into chromosomes, the other(s) can be in the scaffold state). You should know ahead of time **which** genomes you are comparing. If you haven't already, download their genomic FASTA files, which are required.

2. Work from the command line

Most work here is done from the command line, by opening the Terminal application. Commands will be shown here and you may be able to get by with no knowledge of the command line, other than knowing how to start it (by starting the Terminal app on the Mac). However, ideally you should be familiar with some basic command line commands such as *cd*, *ls*, *cp*, *rm*, *mv*, *sudo*, *less*, *head* and *tail*, *man*, and *mkdir*. If you do not know how to use the command line or don't know these basic commands, you may want to take some time to learn a bit. Here are two short tutorial web sites:

<http://www.davidbaumgold.com/tutorials/command-line>

<http://mac.appstorm.net/how-to/utilities-how-to/how-to-use-terminal-the-basics>

Here is a longer tutorial from UC Davis:

http://korflab.ucdavis.edu/Unix_and_Perl/current.html

3. Work in the IGGPIPE main directory unless otherwise instructed

While working from the command line to install IGGPIPE, most of the time you will be in the IGGPIPE main directory, unless instructed otherwise. If you unzipped the IGGPIPE zip file in your Documents folder, you would change into the IGGPIPE directory with this command:

```
cd ~/Documents/IGGPIPE
```

4. Know how to use a plain text editor and have one available

You must have a plain text editor you know how to use. If nothing else, the Mac "TextEdit" program will work (use Plain Text format). The free open-source TextWrangler program is strongly recommended, available from the Apple App Store (Applications, App Store) or from:

<http://www.barebones.com/products/textwrangler>

4. Running IGGPIPE

1. Assign a single capital letter name to each genome

IGGPIPE makes frequent use of a single capital letter to refer to a genome. For example, filenames and tab-separated file column names use such letters. Choose a single capital letter you will use to represent each of your genomes. For example, I used H=Heinz and P=pennellii for some of my testing.

2. Copy allParameters.template file to new file allParameters.XY

IGGPIPE has a number of parameters that must be set to the values you desire. These are contained in plain text files whose name starts with "allParameters". IGGPIPE provides a template file named **allParameters.template** that you copy and edit to create your own "allParameters" files containing the parameter settings you want.

You should make a new parameter file for each different set of settings you want to run with IGGPIPE. For example, each new set of genomes would have a different parameter file. Also, if

RUNNING IGGPIPE TO GENERATE IGG MARKERS

you decided to make two different sets of markers from the same genomes, using two different parameter settings, you might make two different parameter files. Here, we assume to begin with that you are only running your genomes with two genomes and a single set of parameters, and we will name the parameter file `allParameters.XY`, where X and Y are capital letters you chose to represent the two genomes. It is convenient to include those letters in the file name, to help you keep multiple parameter file straight.

Copy `allParameters.template` to `allParameters.XY` (substituting your genome letters for X and Y), using either the Mac's Finder or via the command line. For example:

```
cd ~/Documents/IGGPIPE      (or whatever is appropriate for your system)
cp allParameters.template allParameters.HP    (here, X=H, Y=P)
```

Besides the sample parameter file "`allParameters.template`", there are several more sample "`allParameters`" parameter files in the subfolder *allParameters*. These were used for testing IGGPIPE with various genomes. The subfolder *allParameters* was made to keep the files better organized and keep the root folder less messy. You may want to adopt the same strategy and put your "`allParameters`" files in the *allParameters* subfolder or in a subfolder you make yourself.

3. Open the `allParameters.XY` file in your plain text editor

Open the new **`allParameters.XY`** file created above in your plain text editor for editing. For example, you might open the *nano* text editor from the command line:

```
nano allParameters.XY
```

If you are in a hurry, you don't need to read anything in the file, but can simply **search for** "`#@@`", which are comment lines marking items that may need to be changed.

Each "`#@@`" comment says whether it must be changed, might need to be changed, or probably will never need to be changed, etc. Many of these items typically will never need to be changed, so the actual number of changes that need to be made is smaller than it might first appear. Some of the parameters have already been set correctly during installation of IGGPIPE. However, it is recommended that rather than hurrying, you take time to read through the file, as the comments explain the purpose of each parameter, and you will want to know this information, at least for key parameters, to select the right parameter values for your needs.

4. Search for "`#@@`" and set parameters to desired values

Search for "`#@@`" in the `allParameters.XY` file and check each one to see if it needs to be changed. If so, set it to the value you desire. Parameters you will definitely want to review and consider are:

- a. K
- b. N_GENOMES
- c. GENOME_NUMBERS
- d. GENOME_1, GENOME_2, etc.
- e. LMIN
- f. DMAX

- g. AMIN and AMAX
- h. ADMIN and ADMAX
- i. NDAMIN
- j. OVERLAP_REMOVAL
- k. EPCR_MAX_DEV
- l. EPCR_MAX_MISMATCH and EPCR_MAX_GAPS

After finishing changes, save the modified allParameters.XY file.

5. Check Primer3 settings in primer3settings.txt if desired

The file **primer3settings.txt** contains parameter settings for Primer3, which is used to generate the actual primers. It is possible that you might want to use different Primer3 settings from the defaults listed in this file. If so, edit the file and make the desired changes. For example, you might change the parameters that determine the acceptable *range of primer Tm values*. If you have several different setting values you use, you will probably want to keep a directory of different primer3settings.txt files and copy the needed one to "primer3settings.txt" prior to each run of IGGPIPE. When testing IGGPIPE using the allParameters.test file that was created during installation, always copy primer3settings.default.txt to primer3settings.txt first.

The Primer3 user manual (http://primer3.sourceforge.net/primer3_manual.htm) describes all the Primer3 settings file parameters.

6. Understand use of make and "Makefile" for running IGGPIPE

The IGGPIPE software consists of multiple software applications that progressively analyze the genome sequence data and eventually produce candidate IGG marker primers. The task of running all this software has been automated using a "Makefile", which is a file with that name containing commands formatted correctly for reading the allParameters.XY parameter file and running the software applications. The Makefile is applied by using the application named *make*, which was installed when IGGPIPE was installed, if it didn't already exist.

A big advantage of using "Makefile" and *make* is that if something goes wrong (and by Murphy's law, it probably will), the portion of the work successfully completed is not lost, and does not need to be repeated. This is important because it can take quite a long time to run genomes all the way through the IGGPIPE software. Depending on your computer speed and memory, it can take hours or even days. If an error occurs, *make* will stop, and an error message should be visible on the terminal. After fixing the error, all you have to do resume the pipeline commands from the last successful step is re-enter the same *make* command. *make* knows which step to start at because it knows all the files to be produced by the pipeline, as they are specified in allParameters.XY, and it checks to see if the files exist, and starts at the pipeline step whose output file does not exist. Manually deleting a file will also cause *make* to run the pipeline starting at the step needed to make that file.

You must finish editing the allParameters.XY file before trying to run the pipeline using *make*. If that file is ready to go, you can start running IGGPIPE using the command *make* from the

RUNNING IGGPIPE TO GENERATE IGG MARKERS

command line, with additional command arguments. The first argument that is required is of the form "PARAMS=<allParameters filename>". For example, if your allParameters file is named "allParameters.XY", then the *make* command starts out as *make PARAMS=allParameters.XY*.

The remaining command arguments for the *make* command tell which part of the pipeline to run. If no additional argument is given or if the argument is *ALL*, the entire pipeline is run (or as much of it as is needed to resume where a previous error had halted). However, since the choice of some of the parameters, especially the value of K, can have a strong influence on the number of markers found, it is best to run IGGPIPE a few steps at a time and check the output after those steps before proceeding further. The following sections will guide you in this.

Use this command to get a listing of complete usage information for running *make*:

```
make usage
```

That command will use the *less* command to display file *help.txt*. Press the space bar to move through the text, or press *q* to exit from the help text.

If at any point you want to remove ALL files already generated and start anew, you can do that with this command:

```
make PARAMS=allParameters.XY CLEAN=1 ALL      (replacing with your allParameters name)
```

Running IGGPIPE with a *make* command will usually produce a lot of output on the terminal, and some of this output may be important to examine, especially if an error occurs. Since the output might scroll off the screen and be unavailable, it is a good idea to save it, and this can be done by using the *tee* command along with the *make* command. The *tee* command can write everything that is displayed on the terminal to a file also. You might want to make a folder to contain these "log" files:

```
mkdir logFiles
```

To use *tee*, choose a log file name, let's say *make_kmerIsect_HP14.txt*, and then add at the end of your *make* command line the extra commands *| tee logFiles/make_kmerIsect_HP14.txt*, as in this example:

```
make PARAMS=allParameters.XY kmerIsect | tee logFiles/make_kmerIsect_HP14.txt
```

Then, after *make* finishes, you can examine that log file at any time to see what the pipeline output was, for example:

```
less logFiles/make_kmerIsect_HP14.txt
```

You should use the *tee* command each time you run the pipeline unless you are sure you won't want to reexamine the output later. We will not show the *tee* command in the instructions below, however. It is up to you to decide whether to use it.

We have run IGGPIPE on several different genomes to try to anticipate unusual problems and handle them without error, but there are probably many situations that we haven't yet encountered. If you email us with information about errors and their resolution if you were able to resolve them, we'll try to make improvements to IGGPIPE in error handling and in its input data format flexibility to help future users that encounter the error.

RUNNING IGGPIPE TO GENERATE IGG MARKERS

7. Run IGGPIPE with the command "make PARAMS=allParameters.XY kmerIsect"

The first several steps in the pipeline extract unique k-mers from the FASTA files of the genomes, and intersect these to produce a list of common unique k-mers. To run these steps, use this command:

```
cd ~/Documents/IGGPIPE      (or whatever is appropriate for your system)
make PARAMS=allParameters.XY kmerIsect  (replacing with your allParameters name)
```

or, better yet, log to a file also:

```
make PARAMS=allParameters.XY kmerIsect | tee logFiles/make_kmerIsect_HP14.txt
```

If it completes successfully, the end of the command output will show the message:

```
kmerIsect files are up to date.
```

If it says something else, indicating an error occurred, examine the output carefully and try to diagnose and fix the error, then enter the above *make* command again to retry the failed step.

Assuming the first steps completed successfully, count the number of k-mers in the common unique k-mer file, which is located in the *Kmers* subfolder of the output folder you specified in your *allParameters* file for parameter *DIR_IGGPIPE_OUT*. The name of the file is *isect.kmers* and it is a text file containing one k-mer per line and nothing else. You can look at it with the *less* command if you want to. To count the number of k-mers in it, use the *wc -l* command, which counts lines in a file, like this (replacing *outFolderForMyProject* with your output folder name):

```
wc -l outFolderForMyProject/Kmers/isect.kmers
```

If it shows that you have, say, several million or more, that is good. A few million or less might be too few to generate enough markers. Tens of millions might be too many and cause subsequent pipeline steps to take a very long time. The number of k-mers is influenced by both the value of *K* in the parameter file and by how different the genomes are. Very similar genomes might never have several millions of k-mers. The number of common unique k-mers increases as *K* is increased, up to a point, then starts to decrease. For similar genomes, you might need to edit the parameter file and decrease *K* by 1. With the tomato/pennellii genomes, *K*=14 was a good number. With *Arabidopsis thaliana* Col-0 and Ler-0 accessions, which were much more similar to one another than tomato/pennellii, a value of *K*=13 worked better. For very large genomes, you might need to increase *K* by 1. Since the number of k-mers goes up dramatically with increasing *K*, you probably will never raise *K* above 15, perhaps 16. The default setting in the parameter file for the output directory parameter *DIR_IGGPIPE_OUT* is to include the value of *K* in the directory name. This means you can run IGGPIPE with one value of *K*, then change it and run it again and the output will go into a new directory.

If you feel you have too few or too many k-mers, then you should increase or decrease *K* by 1 and try again. Try both an increase and a decrease in *K* to see how the number of k-mers is affected. Set the parameter file to the value of *K* that seems best to you for proceeding with additional pipeline steps.

8. Run IGGPIPE with the command "make PARAMS=allParameters.XY findLCRs"

The next few steps of the pipeline analyze the common unique k-mers to find locally conserved regions (LCRs). To run these steps, use this command:

RUNNING IGGPIPE TO GENERATE IGG MARKERS

```
make PARAMS=allParameters.XY findLCRs      (replacing with your allParameters name)
```

If it completes successfully, the end of the command output will show the message:

```
findLCRs files are up to date.
```

If it says something else, indicating an error occurred, examine the output carefully and try to diagnose and fix the error, then enter the above *make* command again to retry the failed step.

Assuming the steps completed successfully, count the number of LCRs in the LCRs output file, which is located in the main output directory (set with the DIR_IGGPIPE_OUT parameter). Its name starts with *LCRs_*, and with many command line interfaces you don't need to enter the full name in a command, but can instead enter *LCRs_**. The LCRs file has one LCR per line. You can look at the first few lines with the *head* command or load the file into Excel or a text editor to examine it, if you want. The command line to count the LCRs would look like this:

```
wc -l outFolderForMyProject/LCRs*
```

If you had too few common k-mers you might also have too few LCRs. A million or more LCRs would be nice. The fewer you have, the fewer markers you are likely to get. If there seem to be too few, check the pipeline output. It will show the number of common unique k-mers it processes (it processes them in batches), and the number remaining after it enforces DMIN, LMIN, and KMIN on the reference genome. If these numbers fall dramatically towards 0, it indicates that either there are no good LCRs between the two genomes, or the parameters DMIN, LMIN, and/or KMIN might be too strict. However, expect a pretty big drop with the LMIN step, because typically a large fraction of the common unique k-mers are too close together, with too much separation from the next k-mer, to form a useful LCR.

If you feel you have too few LCRs, then you should try editing the parameter file and changing the DMIN, LMIN, and/or KMIN parameters and try again. The default value for the LCRs_ filename, set by the parameters SFX_LCR_FILE and PATH_LCR_FILE, includes the values of DMIN, LMIN, and KMIN in the filename, so if you change the values, you can simply re-run the pipeline with the same *make* command, and it will generate a new LCRs_ file with a different name, without repeating preceding pipeline steps that do not need to be repeated.

Set the parameter file to the values for DMIN, LMIN, and KMIN that seem best to you and re-run this pipeline step before proceeding with additional pipeline steps.

The columns in the LCRs* file are described in Table 4.

9. Run IGGPIPE with the command "**make PARAMS=allParameters.XY findIndelGroups**"

The next step of the pipeline analyzes the LCRs to find Indel Groups that satisfy the parameters AMIN, AMAX, ADMIN, ADMAX, NDAMIN, and MINFLANK. To run this step, use this command:

```
make PARAMS=allParameters.XY findIndelGroups      (replacing with your allParameters nam
```

If it completes successfully, the end of the command output will show the message:

```
findIndelGroups files are up to date.
```

RUNNING IGGPIPE TO GENERATE IGG MARKERS

If it says something else, indicating an error occurred, then as usual, examine the output carefully and try to diagnose and fix the error, then enter the above *make* command again to retry the failed step.

Assuming the step completed successfully, count the number of Indel Groups in the two output files. One output file includes all Indel Groups found, even when they overlap one another. The other output file includes only non-overlapping Indel Groups, which were determined based on the setting of the parameter `OVERLAP_REMOVAL`. The output files are located in the main output directory and their names start with *IndelGroups* followed by *Overlapping_* and *Nonoverlapping_*. You can look at the first few lines with the *head* command or load the files into Excel or a text editor to examine them, if you want. The command line to count the Indel Groups in both files would look like this:

```
wc -l outFolderForMyProject/IndelGroups*
```

The number of Indel Groups is of the same order of magnitude as the number of markers you will obtain, so pay close attention to the count. If you had too few LCRs you might also have too few markers. Tens of thousands of overlapping markers and thousands of non-overlapping ones are nice numbers to have. If there seem to be too few, you may want to experiment with different values for the `AMIN`, `AMAX`, `ADMIN`, `ADMAX`, `NDAMIN`, and `MINFLANK` parameters. As with the LCRs, these values are normally part of the IndelGroups output file filenames, so if you edit the parameter file and change the values, then re-run the pipeline by entering the same *make* command above, new files are produced with new names.

Set the parameter file to the values for the parameters that seem best to you and re- run this pipeline step before proceeding with remaining pipeline steps.

The columns in the IndelGroups* files are described in Table 3.

10.Run IGGPIPE with the command "**make PARAMS=allParameters.XY findPrimers**"

The next several steps of the pipeline extract DNA sequences for each Indel Group from all genomes using the parameter `EXTENSION_LEN`, then Primer3 (actually, primer3_core) is run to design primers for each Indel Group, using the primer design parameters in file primer3settings.txt.

To run these steps, use this command:

```
make PARAMS=allParameters.XY findPrimers      (replacing with your allParameters name)
```

Depending on the number of Indel Groups, this can take a *long* time. We have seen it take two days to finish running all Indel Groups through Primer3. This is where extra memory can be helpful.

If it completes successfully, the end of the command output will show the message:

```
findPrimers files are up to date.
```

If it says something else, try to diagnose and fix the error as usual, then enter the same *make* command again to retry the failed step.

Assuming the steps completed successfully, count the number of IGG marker primer pairs in the output file, which is located in the main output directory and has a name that starts with *NonvalidatedMarkers_*. You can look at the first few lines with the *head* command or load the files

RUNNING IGGPIPE TO GENERATE IGG MARKERS

into Excel or a text editor to examine them, if you want. The command line to count the primer pairs in the file would look like this:

```
wc -l outFolderForMyProject/NonvalidatedMarkers*
```

Each line of the file contains one pair of primers. Each pair is a candidate IGG marker, but they have not yet been validated using e-PCR, which will mark a few of them as bad and remove them.

The columns in the NonvalidatedMarkers* file are described in Table 1.

11. Run IGGPIPE with the command "make PARAMS=allParameters.XY removeBadMarkers"

The next several steps of the pipeline run e-PCR on each primer pair in the NonvalidatedMarkers* file. All markers are run through e-PCR once for each genome. After that, an R script is run which examines the e-PCR results and removes from the NonvalidatedMarkers* primer pairs all those pairs that failed the e-PCR test in or more genomes, and writes new files with the validated IGG marker primer pairs.

To run these steps, use this command:

```
make PARAMS=allParameters.XY removeBadMarkers      (replacing with your allParameters na
```

Depending on the number of primer pairs, this can take a *long* time. We have seen it take two days *per genome* to finish running all primer pairs through e-PCR. This is where extra memory can be helpful.

If it completes successfully, the end of the command output will show the message:

```
removeBadMarkers files are up to date.
```

If it says something else, try to diagnose and fix the error as usual, then enter the same *make* command again to retry the failed step.

Assuming the steps completed successfully, count the number of validate IGG primer pairs in the two output files. One output file includes all validated primer pairs, even when their amplicons overlap one another. The other output file includes only primer pairs that produce non-overlapping amplicons, determined based on the setting of the parameter OVERLAP_REMOVAL. The output files are located in the main output directory and their names start with *Markers* followed by *Overlapping_* and *Nonoverlapping_*. You can look at the first few lines with the *head* command or load the files into Excel or a text editor to examine them, if you want. The command line to count the primer pairs in both files would look like this:

```
wc -l outFolderForMyProject/Markers*
```

Each line of the file contains one pair of primers. Each pair is an e-PCR-validated IGG marker, essentially the final output of the pipeline.

The columns in the Markers* file are described in Table 1.

12. Run IGGPIPE with the command "make PARAMS=allParameters.XY ALL"

Now you can run the rest of the pipeline to finish up. The final step produces some plots of marker statistics and density. This is done with this command:

RUNNING IGGPIPE TO GENERATE IGG MARKERS

```
make PARAMS=allParameters.XY plotMarkers      (replacing with your allParameters name)
```

However, you can also use this command, which is the command that runs the entire pipeline, all steps from start to finish:

```
make PARAMS=allParameters.XY ALL      (replacing with your allParameters name)
```

This same command could have been used from the very start, to simply run the entire pipeline. It can be used *at any time* to re-run the pipeline. For each step of the pipeline that was previously run successfully, it will simply output a message saying that the files from that pipeline step are up to date. However, for any step that was either not run at all, or failed, it will attempt to re-run that step and following steps that depend on it. At the very end, when it has successfully completed all pipeline steps, it issues this message:

```
ALL files are up to date
```

This is the message you want to see to know that the pipeline has completed successfully.

If you edit the parameter file and change parameters (so that the names of the files that are produced are also changed, since the parameter values are contained within the file names), then when you use the above *make* command, the pipeline re-runs starting at whatever step uses the changed parameters, so you can easily make parameter changes and try again.

Sometimes you may want to force the pipeline to re-run starting at a certain step. For example, maybe you want to re-run the pipeline starting at the *findLCRs* step. To do this, you can either delete the output files produced by that step, or use this command to delete them:

```
make PARAMS=allParameters.XY CLEAN=1 findLCRs      (replacing 'findLCRs' with whatever s
```

After that, you can re-run the pipeline with the *make ... ALL* command shown above to generate the files anew.

The final step, *plotMarkers*, produces several output files in the main output directory. One file is a .pdf file whose name starts with *MarkerCounts_*, containing plots of counts of markers on each chromosome of each genome. You can examine that file with any .pdf file viewer. The other files are .png image files whose names start with *MarkerDensity_* and end with *_X.plot.png*, where X is replaced with the genome letters you assigned for your analysis. Each .png file has an image of the chromosomes with lines showing the positions of each marker. When multiple scaffolds are used, the file limits the output to the first several scaffolds.

13. Comments about pipeline result files

Here is a summary of the filenames produced by each *make* step, in pipeline order:

<i>make</i> command or other command	Produces output file(s)
a. <i>make</i> PARAMS=myFilename getSeqInfo	GenomeData/*.idlens
b. <i>make</i> PARAMS=myFilename getContigFile	GenomeData/*.contigs
c. <i>make</i> PARAMS=myFilename getKmers	Kmers/Kmers_*.kmers
d. <i>make</i> PARAMS=myFilename kmerStats	Kmers/Kmers_*.stats
e. <i>make</i> PARAMS=myFilename sortKmers	Kmers/Kmers_*.sorted

RUNNING IGGPIPE TO
GENERATE IGG MARKERS

<i>make</i> command or other command	Produces output file(s)
f. make PARAMS=myFilename kmerIsect	Kmers/sect.kmers
g. make PARAMS=myFilename getGenomicPos	Kmers/Kmers_*.isect
h. make PARAMS=myFilename mergeKmers	Kmers/Kmers_*.merge
i. make PARAMS=myFilename getCommonUniqueKmers	Kmers/common.unique.kmers
j. make PARAMS=myFilename findLCRs	<i>LCRs_*.tsv, BadKmers_*.tsv</i>
k. make PARAMS=myFilename findIndelGroups	<i>IndelGroupsOverlapping_*.tsv,</i> <i>IndelGroupsNonoverlapping_*.tsv</i>
l. make PARAMS=myFilename getDNaseqs	<i>DNaseqs_*.dnaseqs</i>
m. make PARAMS=myFilename findPrimers	<i>NonvalidatedMarkers_*.tsv</i>
n. make PARAMS=myFilename ePCRtesting	<i>MarkerErrors_*.tsv</i>
o. make PARAMS=myFilename removeBadMarkers	<i>MarkersOverlapping_*.tsv,</i> <i>MarkersNonoverlapping_*.tsv</i>
p. make PARAMS=myFilename plotMarkers	<i>MarkerCounts_*.pdf, MarkerDensity_*.png</i>
q. make PARAMS=myFilename IndelsSNPs	<i>MarkersNonoverlapping_*.indels.tsv,</i> <i>MarkersNonoverlapping_*.snps.tsv</i>
r. make PARAMS=myFilename plotIndels	<i>MarkersNonoverlapping_*.indels.pdf</i>
s. Rscript code/R/dotplot.R dotplot.template	<i>LCRs_*.dotplot.png</i>
t. Rscript code/R/annotateFile.R annotate.template	<i>MarkersAnnotated_*.tsv</i>
u. Rscript code/R/annotateFile.R annotate/ HP11_isInNearColumn.markers	<i>MarkersAnnotated_WithInNearFeatures_*.indels.tsv</i>
v. Rscript code/R/annotateFile.R annotate/ HP11_to_gff3.markers	<i>MarkersAnnotated_GFF3_*.gff3</i>

(Note that some of the files listed above are produced by steps to be described below).

The marker file names and other file names of files in the main output directory are very long and cumbersome, because they include parameter values in them. You may want to copy files to a shorter name to work with them.

The meaning of *overlapping* and *non-overlapping* should be clear from the explanation of the parameter OVERLAP_REMOVAL in the comments in allParameters.XY.

The various .tsv files can be loaded into Excel to examine, and they can also be post-processed (see below) to change them into other formats.

Tables describing each column in each .tsv file type are at the end of this document.

5. Post-processing tools

1. Finding Indels and SNPs

An R program that is NOT run as part of the pipeline when the *make ... ALL* target is built, but which can be run using *make ... IndelsSNPs*, is able to read a file of LCRs, non-overlapping Indel Groups, or non-overlapping Markers, extract the DNA sequences from the genomes in each LCR or Marker region and align them, then locate all Indels and SNPs in the aligned sequences and write their positions to files. The program is called *alignAndGetIndelsSNPs.R*. You should already have set its input file name in your "allParameters" parameter file (PATH_INDELS_SNPS_INPUT_FILE). Run it as follows:

```
make PARAMS=allParameters.XY IndelsSNPs      (replacing with your allParameters name)
```

This produces two files in your output folder whose names end with "indels.tsv" and ".snps.tsv", containing tables of all Indels and SNPs found in the LCR or marker regions. Examine them to see the data they contain. The columns are described in Table 5 and Table 6.

2. Plotting Indel information

Another R program that is NOT run as part of the pipeline when the *make ... ALL* target is built, but which can be run using *make ... plotIndels*, reads the Indels file produced by *make ... IndelsSNPs* and plots information from it in a pdf file. The program is called *plotIndels.R*. Run it as follows:

```
make PARAMS=allParameters.XY plotIndels      (replacing with your allParameters name)
```

This produces a file in your output folder whose name ends with "indels.pdf", containing plots of various Indel information. Examine it to see the plots it contains.

3. Dot plots

The LCRs_*.tsv output file contains locally conserved regions associated with common unique k-mers. It represents a whole genome alignment between the genomes used in IGGPIPE analysis. An R program, *dotplot.R*, is provided that can plot this data as a dot plot.

This program is run by first copying the text file "dotplot.template" to a new name (e.g. *dotplot.XY*) and editing it to specify the parameters of the dot plot. Comments in the file describe each parameter. The program is then run from the command line with a command like this:

```
cd ~/Documents/IGGPIPE      (or whatever is appropriate for your system)
Rscript code/R/dotplot.R dotplot.XY      (or whatever name you gave the parameter file)
```

When it finishes running, the dot plot output file can be found in the place and under the name specified in the parameter file. Use multiple parameter files with different settings to explore different regions of the genomes in greater resolution (parameters include what region of the genome is to be plotted).

The "dotplot.template" file is configured for generating a dot plot file using the LCRs generated via the *allParameters.test* configuration file.

Besides the sample parameter file "dotplot.template" (which has settings for testing the IGGPIPE installation), there are several more sample "dotplot" parameter files in the subfolder *dotplot*, that were used for plotting data from various genomes that IGGPIPE was tested with, although the parameter file is straightforward and you probably don't need other examples to work from. You may want to put your own "dotplot" parameter files in subfolder *dotplot* or your own subfolder to keep them organized.

4. Annotating marker files with other position data and producing GFF3 and GTF files

You may want to make your marker data more conveniently available. For example, you might want to convert it to GFF3 file format so you can add a "marker" track to a genome browser. Or, you may have other genome position data that you would like to have associated with your marker data, such as a file giving positions of introgressions of one genome within another (you might want a column in the marker file showing which introgressions the marker was near). As another example, you might want to add a column in the marker file containing the names of the genes closest to the marker, and the distance to the genes. All of these situations and more can be handled by an R program, *annotateMarkers.R*, provided with IGGPIPE. The program can read and write files of type .tsv (tab-separated variable), .csv (comma-separated variable), .gff3 (general feature format), or .gtf (gene transfer format), all common formats used to hold genome browser track data or FASTA file annotation data. It can add, remove, edit, and rename columns. It can read two separate files and merge their data. It can convert from one of these file formats to another.

This program is run by first copying the text file "annotate.template" to a new name (e.g. *annotateIntrogressions.XY* or *addGeneInfo.XY* or *makeGFF3.XY*) and then editing it to specify the parameters for the annotation and/or file conversion. Comments in the file describe each parameter. The program is then run from the command line with a command like this:

```
cd ~/Documents/IGGPIPE      (or whatever is appropriate for your system)
Rscript code/R/annotate.R addGenes.XY    (or whatever name you gave the parameter file)
```

When it finishes running, the output files can be found in the place(s) and under the name(s) specified in the parameter file.

Besides the sample parameter file "annotate.template" (which has settings for testing the IGGPIPE installation), there are several more sample "annotate" parameter files in the subfolder *annotate*, with file names hinting at what they do, and comments at the start of each file describing what it does. It may be easier to copy one of these and modify it for your needs. You may want to put your own "annotate" parameter files in subfolder *annotate* or your own subfolder to keep them organized.

So, the idea is to use multiple parameter files with different settings to do different types of annotation and file conversion.

Some of the sample parameter files generate .gff3 files that can be added as a track to a genome browser, to display markers in the browser. Instructions for adding the track are given in comments at the start of the parameter file. Two marker files, one for *Arabidopsis thaliana* Col-0 vs. Ler-0 accessions, and the other for *Solanum lycopersicum* vs. *Solanum pennellii* genomes, were created to test IGGPIPE, and the marker files were converted to .gff3 files suitable for making a browser track. These files can be found in subfolders of the *annotate* folder.

File formats can be finicky, especially .gff3 files. An incorrectly formatted file will cause problems with annotateFile.R. When you have problems, if you can submit an issue to the GitHub repository named "BradyLab/IGGPIPE", and attach or insert a copy of your parameter file, that would be helpful. A copy of the input data files would probably also be needed to debug problems, but GitHub does not allow files to be attached. You can email them to us, or find some other way to send them.

5.1. For problems and help:

- Post an issue on GitHub under BradyLab/IGGPIPE repository
- Contact: Ted Toal, twtoal@ucdavis.edu [mailto:twtoal@ucdavis.edu]

6. Tables

Table 1. Columns in MarkersOverlapping_, MarkersNonoverlapping_, NonvalidatedMarkers_, MarkerErrors_ files; X,Y=chosen genome letters

Column	Description
NDA	Number of distinct amplicon sizes, in range NDAMIN..N_GENOMES
Xid	Genome X sequence ID
Xpct	Genome X percent of sequence ID length at which marker is located
XampLen	Genome X amplicon length
Yid	Genome Y sequence ID
Ypct	Genome Y percent of sequence ID length at which marker is located
YampLen	Genome Y amplicon length
YXdif	Difference in length between genomes X and Y amplicons, negative if genome X longer than genome Y
YXphase	Phase of amplicons between genomes X and Y, "+" if both amplicons run in same direction, "-" if opposite directions
prmSeqL	Left side or upstream primer sequence
prmSeqR	Right side or downstream primer sequence
prmTmL	Left side primer Tm
prmTmR	Right side primer Tm
prmLenL	Left side primer length
prmLenR	Right side primer length
XampPos1	Genome X amplicon starting (upstream) position
XampPos2	Genome X amplicon ending (downstream) position, XampPos2 always > XampPos1
YampPos1	Genome Y amplicon starting (upstream) position
YampPos2	Genome Y amplicon ending (downstream) position, YampPos2 > YampPos1 if YXphase is "+", < if "-"

RUNNING IGGPIPE TO
GENERATE IGG MARKERS

Column	Description
kmer1	Common unique k-mer for left side primer region, canonical (lexically smaller of k-mer and its reverse complement)
kmer1strand	N_GENOMES "+" and "-" characters for genomes 1..N_GENOMES. A "+" means k-mer 1 lies on the "+" strand in that genome, "-" means "-" strand.
kmer1offset	Offset in bp of outside (away from amplicon) edge of k-mer 1 from that end of the amplicon. A value of 0 means the amplicon and k-mer ends correspond, >0 means k-mer starts inside the amplicon, <0 means k-mers starts outside it.
kmer2	Common unique k-mer for right side primer region, canonical (lexically smaller of k-mer and its reverse complement)
kmer2strand	Like kmer1strand, for k-mer 2.
kmer2offset	Like kmer1offset, for k-mer 2.
Xseq1	Genome X DNA sequence around left side primer region
Xseq2	Genome X DNA sequence around right side primer region
Yseq1	Genome Y DNA sequence around left side primer region
Yseq2	Genome Y DNA sequence around right side primer region

Table 2. Column reasonDiscarded in MarkerErrors_ files (see Table 1 for other columns)

reasonDiscarded	Description
found multiple	ePCR found multiple amplicons (expected reason)
not found	ePCR didn't find amplicon (should never happen)
wrong seq id	ePCR sequence ID output is wrong (should never happen)
wrong pos	ePCR left and right position output is wrong (should never happen)
wrong posL	ePCR left position output is wrong (should never happen)
wrong posR	ePCR right position output is wrong (should never happen)

Table 3. Columns in IndelGroupsOverlapping_ and IndelGroupsNonoverlapping_ files; X,Y=chosen genome letters

Column	Description
kmer1	Common unique k-mer for left side primer region, canonical (lexically smaller of k-mer and its reverse complement)
kmer2	Common unique k-mer for right side primer region, canonical (lexically smaller of k-mer and its reverse complement)
NDA	Number of distinct amplicon sizes, in range NDAMIN..N_GENOMES
Xid	Genome X sequence ID
Xpos1	Genome X position of upstream end of k-mer 1 on "+" strand
Xpos2	Genome X position of upstream end of k-mer 2 on "+" strand, Xpos1 < Xpos2 always
Xs1	Genome X k-mer 1 strand, "+" or "-"

RUNNING IGGPIPE TO
GENERATE IGG MARKERS

Column	Description
Xs2	Genome X k-mer 2 strand, "+" or "-"
Xctg1	Genome X contig number within sequence Xid of contig containing k-mer 1
Xctg2	Likewise for k-mer 2, Xctg1 = Xctg2 always
XkkLen	Genome X distance from 5' end of k-mer 1 on "+" strand to 5' end of k-mer 1 on "-" strand
Xpct	Genome X percent of sequence ID length at which marker is located
Yid	Genome Y sequence ID
Ypos1	Genome Y position of upstream end of k-mer 1 on "+" strand
Ypos2	Genome Y position of upstream end of k-mer 2 on "+" strand, Ypos1 < Ypos2 if amplicon in X and Y genomes run in the same direction, > if opposite directions
Ys1	Genome Y k-mer 1 strand, "+" or "-"
Ys2	Genome Y k-mer 2 strand, "+" or "-"
Yctg1	Genome Y contig number within sequence Yid of contig containing k-mer 1
Yctg2	Likewise for k-mer 2, Yctg1 = Yctg2 always
YkkLen	Genome Y distance from 5' end of k-mer 1 on "+" strand to 5' end of k-mer 1 on "-" strand
Ypct	Genome Y percent of sequence ID length at which marker is located

Table 4. Columns in LCRs_ and BadKmers_ files; X,Y=chosen genome letters

Column	Description
(none, row name)	Common unique k-mer, canonical representation (the lexically smaller of k-mer and its reverse complement)
X.seqID	Genome X sequence ID
X.pos	Genome X position of upstream end of k-mer on "+" strand relative to start of X.seqID
X.strand	Genome X k-mer strand, "+" or "-"
X.contig	Genome X contig number within sequence X.seqID sequence of contig containing the k-mer
X.contigPos	Genome X position of upstream end of k-mer on "+" strand relative to start of X.contig
Y.seqID	Genome Y sequence ID
Y.pos	Genome Y position of upstream end of k-mer on "+" strand relative to start of Y.seqID
Y.strand	Genome Y k-mer strand, "+" or "-"
Y.contig	Genome Y contig number within sequence X.seqID sequence of contig containing the k-mer
Y.contigPos	Genome Y position of upstream end of k-mer on "+" strand relative to start of Y.contig
LCR	Integer LCR number to which this k-mer is assigned, each LCR has a unique LCR number assigned to it

Table 5. Columns in *.indels.tsv files; X,Y=chosen genome letters

Column	Description
ID	Unique ID tying row back to originating input file row. LCR input files: LCRnumber. IndelGroup and Markers files: refID_refPos1_refPos2.
phases	Phase of each genome incl. ref. genome, relative to ref. genome, string of +/- chars, + : same direction, - : opposite direction.",
idx	Starts at 1 and counts each Indel within an ID. For given ID (input row), number of Indels in that region is max idx value. If more than two genomes, entire region where alignment has a gap in one or more genomes is counted as one Indel even if multiple gap regions occur in different genomes.
Xdel,Ydel	Total number of deleted bps within the Indel in genomes X,Y. With 2 genomes, del = 0 in genome with insertion (no gaps), del > 0 in genome with deletion (gaps). With >2 genomes, del can be non-zero for all genomes. A genome has only insertions in the Indel if del is 0, and it has only deletions if end-start-1 = 0, and otherwise it has a mixture of at least one insertion and one deletion within the Indel interval.
Xid,Yid	Sequence ID of the Indel in genomes X,Y.
Xstart,Xend,Ystart,Yend	Overall Indel starting and ending position in genomes X,Y. start/end are positions of bps just BEFORE first and AFTER last Indel gap in any genome, so they refer to the same two bps in all genomes. Always start < end. If - phase, start is bp just AFTER, end is bp just BEFORE, opposite of +. Length of the Indel region in each genome is end-start-1.

Table 6. Columns in *.snps.tsv files; X,Y=chosen genome letters

Column	Description
ID	Unique ID tying row back to originating input file row. LCR input files: LCRnumber. IndelGroup and Markers files: refID_refPos1_refPos2.
phases	Phase of each genome incl. ref. genome, relative to ref. genome, string of +/- chars, + : same direction, - : opposite direction.",
idx	Starts at 1 and counts each SNP within an ID. For given ID (input row), number of SNPs in that region is max idx value.
Xid,Yid	Sequence ID of the SNP in genomes X,Y.
Xpos,Ypos	SNP position in genomes X,Y.
Xval,Yval	SNP value in genomes X,Y.